

# JIN LAB Nuclei extraction protocol (from perinatal mouse brain)

9/1/2023 / Graham Anderson / v1

## NOTES

Optimal douncing and lysis conditions depend on the specific tissue used and on intended downstream processes. Amount of strokes using dounce pestles A or B, post-douncing lysis rest, and formulation of NIB-NP and NIB-BSA may all be optimized for different tissues. Those conditions here are intended for fresh P3-P4 mouse cortex, a fragile specimen.

## MATERIALS

Dounce Homogenizer (Sigma-Aldrich D8938)

70um Cell Strainers (Corning 352350)

**Ultra Pure Water** (Fisher Scientific RGF-3410)

**10X PBS** (Life Technologies AM9625)

**NIM (Nuclei Isolation Medium)** ([fresh](#) or can be stored overnight at 4C)

- Sucrose solution (1.5M, or make fresh from solid) (Sigma-Aldrich 84097)
- KCl (2M) (Life Technologies AM9640G)
- HEPES (1M, pH 7.5) (ThermoFisher 15630080)
- MgCl<sub>2</sub> (1M) (Life Technologies AM9530G)

**NIB-NP (Nuclei Isolation Buffer with NP40)** ([fresh](#))

- NIM
- Dithiothreitol (DTT) (1M) (Sigma-Aldrich 646563)
- Kollidon, or Poly(1-vinylpyrrolidone-co-vinyl acetate) (Sigma-Aldrich 190845)
- Nonidet P40 (NP-40), or substitute/alternative (10% v/v) (MilliporeSigma 492016)

**NIB-BSA (Nuclei Isolation Buffer with BSA)** ([fresh](#))

- NIM
- DTT
- Bovine Serum Albumin (BSA) (solid or 10% w/v) (Sigma-Aldrich B6917)

## METHODS

1. If need to prepare NIM, thaw or prepare **1.5M Sucrose** according to **Table 1**.
  - a. Once dissolved, filter with sterile 0.22um filter. Excess can be aliquoted and stored at -20C. Check for white bacterial growth before using thawed sucrose.

**Table 1: Sucrose Prep (more than 1/2 as much as NIM, filtration will cause loss)**

Final [Sucrose] (M)	Mass Sucrose (g)	Total Volume (mL)
1.5	20.538	40

2. Prepare **NIM** according to **Table 2**. Vortex to mix and then filter with 0.22um filter.

**Table 2: NIM Prep (6mL per sample preparation; more for 1 sample - filtration will cause loss)**

Reagent		Final Concentration (mM)	Amount to Add	
Name	Initial Concentration (M)		#	Units
Sucrose	1.5	250	1	mL
KCl	2	25	75	<b>uL</b>
HEPES	1	10	60	<b>uL</b>
MgCl <sub>2</sub>	1	5	30	<b>uL</b>
Water to fill			4.835	mL
<b>Total</b>			6	mL

3. Thaw or dissolve DTT and BSA.
4. Prepare **NIB-NP** and **NIB-BSA** on ice according to **Tables 3** and **4**, keep chilled. Dissolving Kollidon may take some time.

**Table 3: NIB-NP Prep (4mL per sample preparation)**

Reagent	Initial Concentration		Final Concentration		Amount to Add	
	Amount	Units	Amount	Units	#	Units
NIM					3.956	mL
DTT	1	M	1	mM	4	<b>uL</b>
Kollidon	100	% (solid)	10	mg/mL	40	<b>mg</b>
NP-40	10	% (liquid)	0.1	%	40	<b>uL</b>
<b>Total</b>					4	mL

**Table 4: NIB-BSA Prep (at least 2mL per sample preparation)**

Reagent	Initial Concentration		Final Concentration		Amount to Add	
	Amount	Units	Amount	Units	#	Units
NIM					1798	uL
DTT	1	M	1	mM	2	uL
BSA*	10	% (liquid)	1	%	200	uL
<b>Total</b>					2	mL

\*Prepare 10% by dissolving 0.1g solid filled up to 1mL water

5. Fill dounce vials with NIB-NP to the line (~1.6mL) and embed in ice (past the fill level). Embed centrifuge tubes in ice. Ensure the centrifuge is cooling to 4C. Nuclei should always be kept at 4C/on ice.
6. Transfer tissue to dounce vial and dounce 3x with pestle A, then 3x with pestle B.
  - a. Do not remove vials from ice to dounce. Bring rods only close to bottom of vials or gently touch it, do not add pressure or grind.
  - b. The number of douncing and post-dounce lysis rest time will need to be adjusted for different specimen types (age, tissue types, fresh versus frozen samples, etc).
7. Transfer the lysis to a fresh tube, **add 2 mL fresh** NIB-NP, and continue lysing for 5 minutes resting embedded in ice.
8. Transfer the sample through a pre-wetted, 70 um strainer to a centrifuge tube.
 

Note: pre-wet the strainer with PBS to prevent sample loss.
9. Centrifuge at 200g at 4C for 8 minutes.
 

Note: shorter centrifugation time may lead to loss of nuclei.
10. Carefully remove supernatant without disrupting pellet, then resuspend by thoroughly but gently pipetting in 0.5mL **NIB-BSA**.
11. Repeat centrifugation and resuspension in step 10. Sample is now ready for FACS or other analysis. Check the nuclei morphology under the scope to ensure the integrity of the sample.